



BEAMSTOP	MD2S or MD3 type	Beamstop size	Range: Sample / beamstop
Ref .	MD2S or MD3 standard	D (mm)	A (mm)
300-500-60	300	0,50	60 - 8
300-500-80	300	0,50	80 - 28
300-400-60	300	0,40	60 - 8
300-400-80	300	0,40	80 - 28
300-300-60	300	0,30	60 - 8
300-300-80	300	0,30	80 - 28

MD2S or MD3 type	Ref.
MD2S or MD3 standard	300
MD2S or MD3 mirrored	400



/	/	Contrôle après montage	WC	Présence de colle
Rep	Nb	NOTA	Matière	Observations
C				
B				
A	23/10/15	First Issue		
Ind.	Date	Modification		
Service : Instrum.Scientif		Tol. gén. ISO:	Js9	
Verif : RS		Tol. gén. Rugosité	Ra 1.6	
Dess: CHU		Date: 23/10/2015	A2	
Ech: 3:1		Sous Ens: Organs		Projet: MD2S - MD3
TITRE:				N°plan
Beamstop MD2S - MD3				CME-MD2S-1010-A

